-  **(EN)** Flip-Top Dilution Bottle with Buffered Peptone Water (ISO)
-  **(FR)** Flip-Top avec Eau Peptonée Tamponnée (ISO)
-  **(DE)** Flip-Top Verdünnungsröhrchen mit gepuffertem Peptonwasser (ISO)
-  **(IT)** Flip-Top flacone per diluizione con acqua peptonata tamponata (ISO)
-  **(ES)** Flip-Top Frasco de Dilución con Agua Peptonada Tamponada (ISO)
-  **(NL)** Flip-Top, met gebufferd pepton water (ISO)
-  **(SV)** Flip-Top Dilution Bottle med buffrat peptonvatten (ISO)
-  **(DA)** Flip-Top fortyndingsflaske med Buffered Peptonvand (ISO)
-  **(NO)** Fortynningsflaske med vippelokk bufret peptonvann (ISO)
-  **(FI)** Flip-Top Laimennuspullo puskuroidulla peptonivedellä (ISO)
-  **(PT)** Flip-Top Frasco de Diluição com Água Peptonada Tamponada (ISO)
-  **(EL)** Flip-Top Φιάλη Αραίωσης με Ρυθμιστικό Υδατικό Διάλυμα Πεπτόνης (ISO)
-  **(PL)** Butelka do rozcieńczeń typu Flip-Top ze zbuforowaną wodą peptonową (ISO)
-  **(HU)** Felpattintható-kupakú higító palack pufferolt pepton vízzel (ISO)
-  **(CS)** Ředicí nádobka s uzávěrem Flip Top a pufrovanou peptonovou vodou ISO
-  **(RO)** Sticlă de diluție Flip-Top cu apă peptonată tamponată (ISO)
-  **(RU)** Флакон Флип-Топ с забуференной пептонной водой (ISO)
-  **(TR)** Flip-Top Seyreltme Şişesi Tamponlanmış Peptonlu Su ile (ISO)
-  **(JA)** フリップトップ 希釈水 緩衝ペプトン水 (ISO)
-  **(ZH)** Flip-Top 缓冲蛋白胨水稀释瓶 (ISO)
-  **(TH)** ขวดบรรจุสารละลายบัฟเฟอร์เปปโตนวอเตอร์ (ISO)
-  **(KO)** Flip-Top 희석액 (Buffered Peptone Water (ISO))

Product Instructions

Flip-Top Dilution Bottle with Buffered Peptone Water (ISO)

Description

The 3M™ Flip-Top Dilution Bottle with Buffered Peptone Water (ISO) is a diluent for the enumeration of microorganisms from food, beverage, and environmental samples. The diluent conforms to the formula specified in the International Organization for Standardization (ISO), ISO 6579¹ and is as follows:

Enzymatic digest of casein 10.0 g

Sodium chloride 5.0 g

Disodium hydrogen phosphate dodecahydrate ($\text{Na}_2\text{HPO}_4 \cdot 12\text{H}_2\text{O}$) 9.0 g

Potassium dihydrogen phosphate (KH_2PO_4) 1.5 g

Water 1000 mL

pH is adjusted as necessary, so that after terminal gamma irradiation, it is 7.0 ± 0.2 at 25°C.

3M Food Safety is certified to ISO 9001 for design and manufacturing.

Safety

The user should read, understand, and follow all safety information in the instructions for the 3M Flip-Top Dilution Bottles. Retain the safety instructions for future reference.

NOTICE: Indicates a hazardous situation, which, if not avoided, could result in property damage.

NOTICE

- 3M has not documented 3M Flip-Top Dilution Bottle with Buffered Peptone Water (ISO) diluent for use in industries other than food. For example, 3M has not documented 3M Flip-Top Dilution Bottle with Buffered Peptone Water (ISO) for testing water, pharmaceuticals or cosmetics.
- 3M Flip-Top Dilution Bottle with Buffered Peptone Water (ISO) has not been tested with all possible food products, food processes, testing protocols or with all possible strains of bacteria.
- Do not use 3M Flip-Top Dilution Bottle with Buffered Peptone Water (ISO) in the diagnosis of conditions in humans or animals.
- The user must train its personnel in proper testing techniques: for example, Good Laboratory Practices², ISO 7218³, or ISO 17025⁴.
- Select commercially made buffered peptone water media formulated to meet the requirements of ISO 6887 such as 3M Flip-Top Dilution Bottle with Buffered Peptone Water (ISO), may inhibit the DNase reaction resulting in no pink zone formation when the 3M™ Petrifilm™ Staph Express (STX) Plate is used with the 3M™ Petrifilm™ Staph Express (STX) Disk.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.3M.com/foodsafety, or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

Limitation of 3M Liability

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage

Upon receipt, store unopened 3M Flip-Top Dilution Bottle with Buffered Peptone Water (ISO) at temperatures 20-25°C.

△ Disposal

After use, 3M Flip-Top Dilution Bottle with Buffered Peptone Water (ISO) may contain microorganisms that may be a potential biohazard. Follow current industry standards or local regulations for disposal.

Instructions for Use

If used with 3M™ Petrifilm™ Plates:

Plating

Refer to the appropriate 3M Petrifilm Plate package insert for plating instructions.

Incubation

Refer to the appropriate 3M Petrifilm Plate package insert for Incubation instructions.

Interpretation

Refer to the appropriate 3M Petrifilm Plate package insert for Interpretation instructions.

For further information refer to the appropriate 3M Petrifilm Plate "Interpretation Guide." If you have questions about specific applications or procedures, please visit our website at www.3M.com/foodsafety or contact your local 3M representative or distributor.

References

1. ISO 6579. Microbiology of food and animal feeding stuffs Horizontal method for the detection of *Salmonella* spp.
2. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

Refer to the current version of the standard methods listed above.

Explanation of Symbols

www.3M.com/foodsafety/symbols

3M Food Safety

3M United States

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-563-2921

3M Latin America

3M Center
Bldg. 275-5W-05
St. Paul, MN 55144-1000
USA
1-954-340-8263

3M Europe and MEA

3M Deutschland GmbH
Carl-Schurz-Strasse 1
D41453 Neuss/Germany
+49-2131-14-3000

3M United Kingdom PLC

Morley Street,
Loughborough
Leicestershire
LE11 1EP
United Kingdom
+(44) 1509 611 611

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Euro Plaza
Gebaude J, A-1120 Wien
Kranichberggasse 4
Austria
+(43) 1 86 686-0

3M Asia Pacific

No 1, Yishun Avenue 7
Singapore, 768923
65-64508869

3M Japan

3M Health Care Limited
6-7-29, Kita-Shinagawa
Shinagawa-ku, Tokyo
141-8684 Japan
81-570-011-321

3M Australia

Bldg A, 1 Rivett Road
North Ryde, NSW 2113
Australia
61 1300 363 878

3M Health Care

2510 Conway Ave
St. Paul, MN 55144 USA
www.3M.com/foodsafety

© 2018, 3M. All rights reserved.
3M and Petrifilm are trademarks of 3M. Used under license in Canada.
34-8722-6268-7

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Referenties

1. ISO 6579. Microbiology of food and animal feeding stuffs Horizontal method for the detection of *Salmonella* spp.
2. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
3. ISO 7218. Microbiology of food and animal feeding stuffs - General requirements and guidance for microbiological examinations.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

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1. ISO 6579. Microbiology of food and animal feeding stuffs - Horizontal method for the detection of *Salmonella* spp.
2. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
3. ISO 7218 Microbiology of food and animal feeding stuffs - General requirements and guidance for microbiological examinations.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

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3M Österreich GmbH