

Isolation JPAJASWADS www.mwe.co.uk

Unique self-contained detection for Listeria, Salmonella, Coliform and Hygiene

STERI

SWAB

Ш

STERII 05A20

JAN 06

STERI

04L09

NOV₀₆

TRANSWAB

<

Y Ш

LIST

site

Sampling

Locatio

Each test is entirely self-contained **Rapid detection of pathogens** Easy to use & easy to interpret Organisms remain viable for further testing, confirmation, investigation, etc





Isolation Transwabs®



Coliform Isolation Transwab®

For the early detection of coliforms and E.

coli in the environment.

Escherichia coli and Coliforms are traditionally used as indicator organisms for faecal contamination in water and other environmental samples. Detection of these organisms usually points to poor hygiene at some stage in the production process or pollution of water at source.

Based on a traditional test for Coliforms, the **Coliform Isolation Transwab®** provides a simple method for detecting these organisms - showing a definite colour change from purple to yellow when coliforms are present.

Listeria Isolation Transwab[®]

Easy to use early warning test for Listeria contamination in food production

premises.

Easy-to-use and interpret. Simple colour change indicates contamination. No colour change means it's clear! Self-contained swab test. No other reagents to add. Ideal for control of critical points in the supply and manufacturing chain.





Hygiene Swab™

Same day results

Same day results from Hygiene Test make it possible for food manufacturers to avert potential major production loss by taking appropriate cleaning down procedures to reduce contamination of finished products.

Rapid Screening

Hygiene Test is a rapid test which will detect gross bacterial and fungal contaminations of work surfaces, equipment machinery or other sampling sites.

Helps improve hygiene standards

Regular use of Hygiene Test ensures the maintenance of high standards of general cleanliness in food manufacturing plants and catering establishments.

Order Details		
Product Ref.		
MW503/25		
MW504/25		
MW570/25		

Description	Unit size
Coliform Isolation Transwab [®]	25
Hygiene Swab™	25
Listeria Isolation Transwab $^{^{ extsf{B}}}$	25

Medical Wire & Equipment

 Leafield Industrial Estate,

 Corsham, Wiltshire, SN13 9RT

 Tel: 01225-810361

 Fax: 01225-810153

E-mail: info@mwe.co.uk www.mwe.co.uk

